

CD 2000



Bipolar Discharge System

Model Variations



CD 2000 Type A
Bipolar discharge unit with lower mixed air flow



CD 2000 Type B
Bipolar discharge unit with higher mixed air flow

Description



Fig. 1: CD 2000 The bipolar corona discharge CD 2000 is used in order to neutralize the charge of a dust or liquid aerosol or even to electrically charge it, as required. For this purpose, a current is set on the positive and negative high voltage source so high that it is still fully conducted through the corona discharge generated at the tip in the corona chamber, but no (or few) sparks are generated. This current is programmed via the current limitation of the HV power supplies. The voltage limitation is set slightly higher than is necessary for the current determined, such that this current flows under all circumstances (temperature, condition of the electrode). The positive and the negative high voltage supplies and the CD 2000 are independent of one another and are independently adjustable. **Functional principle**




Fig. 2: Functional principle of corona discharge A dust or liquid aerosol is supplied at the aerosol inlet. The particles in this aerosol have to be discharged. The mixed air inlet is fed with dry pressurized air. The mixed flow air is passed through two ionization chambers. In these chambers the air is enriched with ions from the corona discharge with positive and negative high voltage. The aerosol air stream and the mixed air enriched with ions are mixed in a small mixing chamber. The charges of the particles in the aerosol are neutralized by the ions in the mixed air. This takes a certain time. That means that the complete discharge only takes place in the volume after CD 2000 outlet (e.g in the tubing system or in the upper part of a filter holder). An approximate value for the dwell time is 0.5 seconds. A satisfactory discharging effect was found in configurations that ensure a dwell time of at least 0.2 seconds. **CD 2000 measurement data**




Fig. 3: Fractional separation efficiency with and without CD 2000 With a filter medium or filter (exception: Electret filter medium or Electret filter) the fractional separation efficiency is normally measured higher without a discharge unit than with one. **Attention** Many standards, e.g. EN 779 for room filters, require the use of a discharge unit.

Benefits

- No operation license required for radioactive instruments
- Bipolar discharge through negative and positive ions
- Applicable for solid and liquid aerosols
- Robust design
- Simple operation
- Reliable function
- Low maintenance
- Reduces your operating expenses

Datasheet

Parameter	Description
Power consumption	50 W
Power supply	115 - 230 V, 50 - 60 Hz
Reported data	Voltage: 0 - 6,000 V $\hat{=}$ 0 - 10 V, power: 0 - 1,000 μ A $\hat{=}$ 0 - 10 V
Aerosol outlet connection	Aerosol and fed mixed air, $\varnothing_{\text{inside}} = 12$ mm, $\varnothing_{\text{outside}} = 16$ mm
Special features	Positive and negative high voltages are provided by two independent power supplies, maximum voltage: $\pm 6,000$ V, maximum power: $\pm 1,000$ μ A
Volume flow (suction flow)	0 - 4 m ³ /h
Mixed air connection	Cleaned pressurized air, type A: $\varnothing_{\text{inside}} = 6$ mm, $\varnothing_{\text{outside}} = 8$ mm, type B: $\varnothing_{\text{inside}} = 13$ mm
Operation principle	Ionization with corona
Mains fuse	F 3,15 A, 250 V
Volume flow (mixed air)	Type A: for 2 - 18 m ³ /h, type B: for 3 - 36 m ³ /h
Aerosol inlet connection	$\varnothing_{\text{outside}} = 8$ mm, $\varnothing_{\text{inside}} = 6$ mm

Applications

- Discharge of electrically charged aerosols
- Aerosol research
- Filter testing

Palas GmbH
Partikel- und Lasermesstechnik
Greschbachstrasse 3 b
76229 Karlsruhe
Germany

Managing Partner:
Dr.-Ing. Maximilian Weiß
Commercial Register:
register court: Mannheim
company registration number: HRB 103813
USt-Id: DE143585902



Contact: E-Mail: mail@palas.de Internet: www.palas.de Tel: +49 (0)721 96213-0 Fax: +49 (0)721 96213-33