

MSS 08



The Multi Sampling System MSS 08 enables the connection of central measuring devices (e.g., UF-CPC¹ or ENVI-CPC²) up to eight different sampling ports and easy switching between them. This can be useful, e.g., for nanoparticle concentration monitoring in different rooms or along a process with a single measuring device.

BENEFITS

- Automatic change between 8 measuring positions
- Reduces acquisition and operating costs

APPLICATIONS

- Air quality surveillance in different rooms
- Process monitoring on multiple positions

¹condensation particle counter UF-CPC: <https://www.palاس.de/en//en/product/ufcpc>

²condensation particle counter ENVI-CPC: <https://www.palاس.de/en//en/product/envicpc>

DATASHEET

Volume flow	0 – 10 l/min
Power supply	115 – 230 V, 50/60 Hz
Response time	3 s (for change of position)
Aerosol inlet connection	8 mm