


PDANALYZE



The PDAnalyze software is designed to acquire measurement values and data analysis by all stand-alone particle measuring devices. It enables quick and easy analysis of particle sizes and concentrations using the Promo[®] System, Fidas[®] System, UF-CPC System, U-SMPS System, and U-RANGE System.

BENEFITS

- Optimal information content
 - Display of distributions in the diagram and table form: Number, area, volume, concentration (number and mass)
 - Time sequence of fine dust values
 - Comparison of multiple distributions in one diagram
 - Display, time sequence, and analysis of 24 statistical values
- Acquisition and analysis in as little as 1 s
- Rapid incorporation due to uniform interface, even with various measuring devices
- Color plotting for the representative display of measurement results
- Clear structure
- Regular updates via the Internet free of charge
- Easy data export to text file that can be imported into Excel.
- Easy installation on the PC
- Unlimited copies for various workstations within a company
- Easy to operate
- Reliable function
- Reduces your operating expenses

APPLICATIONS

- Particle measurement and particle size analysis
- Comparing distributions from 5 nm to 100 µm
- Evaluation of emission and immission measurements
- Analysis of long-term measurements over days or months
- Analysis of measurements by optical aerosol spectrometers, fine dust measuring devices, nanoparticle measuring devices
- Optical processing of measurement data for use in presentations

DATASHEET

| | |
|-------------------|--|
| Operating system | Windows® 7, Windows® 10 |
| Processor | Min. Pentium I5, 2000 MHz |
| User memory | From 4 GB |
| Screen resolution | 15.6" (min.1600 • 900 for Notebooks), 24" (min. 1920 • 1080 for PCs) |
| Messgeräte | CPC System, U-SMPS System, Promo® System, Fidas® System |