



The sampling tubes PS 1000 / PS 2000 and PS 3000 can be heated up to 120 °C and easily adapted to the heated welas® aerosol sensors. The heating control is part of the delivery content.

## 优势

- Actively heated, flexible sampling tube
- Easy use with welas® aerosol sensors
- Avoiding of condensation effects during aerosol sampling by heating to the appropriate temperature
- Selection for temperatures up to 120 °C or up to 250 °C

## 应用领域

- Isothermal sample taking of hot aerosols
- Blow-by measurements

## 技术数据

电源	115 V or 250 V
Temperature regulation	Heatable to 120 °C or 250 °C
Diameter	Inside = 6 mm, outside = 8 mm
Dimensions	1 m, 2 m or 3 m (length)