

# PLG 1000



The PLG 1000 is a cold atomizer intended for use in air-conditioned rooms. If the room cannot be air-conditioned, a heatable device version should be used, e.g., PLG 2000 H.

## 优势

- Excellent short-term and long-term dosing constancy
- Best reproducibility with respect to particle size distribution and particle concentration
- Large mass volume range (very low and very high)
- Robust design (optionally resistant against chemically aggressive liquids)
- Compact and light
- Easy to operate, proven in industrial applications
- Reliable function
- Low maintenance

## 应用领域

- Filter industry/oil separators
  - Determination of separation efficiency
  - Determination of fractional separation efficiency
  - Loading test
  - HEPA/ULPA filters
- Test of cooling lubricant separators
- Comparison of particle measurement devices
- Tracer particles
- Flow visualization

## 技术数据

体积流量	1 – 23 l/min
Mass flow (particles)	2.5 g/h
Filling quantity	70 ml
Aerosol outlet connection	$\varnothing_{\text{innen}} = 11 \text{ mm}, \varnothing_{\text{au\ss en}} = 14 \text{ mm}$
Mean particle diameter (number)	0.4 $\mu\text{m}$ (DEHS)
Dimensions	280 • 130 • 100 mm (H • W • D)
重量	Approx. 2 kg